



# DMQ<sup>®</sup>

## damp mop neutral disinfectant cleaner

DMQ is a non-alkaline disinfectant cleaner concentrate specially designed to clean, deodorize and disinfect high gloss floors. Effective as a sanitizer in 30 seconds, DMQ retains the beauty of high gloss floors with the added effectiveness of a disinfectant! Formulated for use in daily maintenance programs, DMQ is a proven one-step cleaner disinfectant in the presence of moderate amounts of organic material.



Providing Quality Maintenance and Chemical Specialty Solutions with Unparalleled Customer Service



Spartan Chemical Company, Inc.  
 1110 Spartan Drive  
 Maumee, OH 43537  
 1-800-537-8990  
 www.spartanchemical.com

**EXCELLENT CLEANING PERFORMANCE**

DMQ is a concentrate containing only the highest quality raw material ingredients. An effective disinfectant when diluted 2 oz. per gallon of water, DMQ is economical enough to use every day. This feature is especially important when DMQ is the disinfectant damp mop choice for ultra high speed floor maintenance programs.

**KILLS THE GERMS... NOT THE SHINE!**

The DMQ formula has a low pH and special non-alkaline ingredients. It will not dull, haze or streak high gloss floor finish. Suitable for the Bloodborne Pathogen Standard.

**GREAT FOR HIGH SPEED FLOOR FINISH PROGRAMS**

DMQ's film-free characteristic plus its germ-killing properties make it especially desirable for the daily damp mopping procedure in ultra high speed floor care programs. Schools, office buildings and other commercial areas where there is a high degree of traffic will particularly appreciate the benefits of a film-free disinfectant damp mop cleaner.

**DMQ IS BACTERICIDAL**

Effective disinfectant against the following bacteria:

- Staphylococcus aureus
- Salmonella enterica

**DMQ KILLS ANTIBIOTIC-RESISTANT BACTERIA\*\***

Effective disinfectant against the following antibiotic-resistant bacteria\*\*:

- Methicillin-Resistant Staphylococcus aureus (MRSA)
- Vancomycin-Resistant Enterococcus faecalis (VRE)

**DMQ KILLS 99% OF BACTERIA\*\*\* IN 30 SECONDS**

Effective sanitizer in 30 seconds against the following bacteria\*\*\*:

- Staphylococcus aureus
- Enterobacter aerogenes

**DMQ IS VIRUCIDAL\***

Effective virucide against the following viruses\*:

- HIV-1 (AIDS Virus)
- Hepatitis B
- Hepatitis C
- Herpes simplex virus Type 2
- Influenza A2/Hong Kong Virus

**AREAS FOR USE**

DMQ is recommended to clean, disinfect, sanitize and deodorize floors, walls, countertops, fiberglass, synthetic marble, stainless steel, glass, vinyl and other hard, non-porous environmental surfaces in:

- |                  |                        |
|------------------|------------------------|
| Schools          | Churches               |
| Public Buildings | Factories              |
| Restaurants      | Manufacturing Sites    |
| Bars             | Food Processing Plants |
| Cafeterias       | Apartment Buildings    |
| Office Buildings | Recreation Facilities  |
| Hotels           | Public Restrooms       |
| Motels           |                        |

**See product label for complete directions for use and more applications.**

**Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting [www.spartanchemical.com](http://www.spartanchemical.com).**

**SPECIFICATIONS**

pH (Concentrate): 5.0-6.0

Color: Yellow

Scent: Light citrus

Stability:

- Shelf @ 24°C/75°F—One year minimum
- Freeze/Thaw—Three cycles

EPA Reg. No. 5741-20



Distributed by:

**PACKAGING**

106265	330-gallon tote
106255	55-gallon drum
106230	30-gallon drum
106215	15-gallon drum
106205	5-gallon pail
106204	1-gallon (4 per case)



Secondary labels are available in English, Spanish, and French.

**GUARANTEE:** Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of the date of manufacture. Use product as directed and read all precautionary statements.

For institutional and industrial use only.

© SCC 2/14 L1062