



Spartan Chemical Company, Inc., 1110 Spartan Drive, Maumee, OH 43537
 Toll Free: 1-800-537-8990, Fax: 419-536-8423, Fax-On-Demand: 888-383-4666
 E-Mail: spartan@spartanchemical.com, Web Site: www.spartanchemical.com
 ©SCC 2003



Iodine & Betadine Remover

PRODUCT DESCRIPTION:

Clothesline Fresh Iodine & Betadine Remover removes iodine, betadine, and chlorine stains from fabrics. Clothesline Fresh Iodine & Betadine Remover will also remove these stains on carpeting! This stain removing power on fabrics and carpet makes Clothesline Fresh Iodine & Betadine Remover ideal for hospitals, nursing homes, and other health care facilities.

Use Clothesline Fresh Iodine & Betadine Remover to remove these spots:

- iodine
- betadine
- chlorine

COLOR AND NUMBER SYSTEM:

Each Clothesline Fresh laundry product has been color coded and number coded for easy product identification. All product information and literature, plus, each product label, carton and secondary label reflects each individual product's color and number, making training easier and faster.

STORAGE:

Clothesline Fresh laundry products can be stored in a wide variety of environments and locations. Clothesline Fresh Iodine & Betadine Remover can withstand three freeze/thaw cycles without affecting the performance of the product.

DIRECTIONS FOR USE:

DO NOT MIX WITH OTHER CHEMICALS.

Some individuals may be sensitive to ingredients in this product. Before use, read product label and MSD sheet.

- Always test fabric for color fastness and compatibility in an inconspicuous area.
 - Before using a specific spotter, attempt to identify the stain, to ensure efficient stain removal.
 - Always wash and rinse fabric thoroughly after spot treatment.
1. Spray Clothesline Fresh Iodine & Betadine Remover liberally onto suspected stains.

2. Allow to remain on stained fabric for:
 - 10 – 15 minutes for iodine stains
 - 2 – 3 minutes for betadine stains
 - 1 – 5 minutes for chlorine stains
3. Wash and rinse fabric thoroughly.

SPECIFICATION DATA:

pH – 9.0 – 9.5

Specific Gravity – 1.022 @ 24°C/75°F

Density – 8.505

Stability -

- a. Shelf – One year minimum @ 24°C/75°F
- b. Accelerated – 30 days minimum @ 49°C/120°F
- c. Freeze/Thaw – Can withstand three freeze/thaw cycles

Phosphate free

Nonylphenol ethoxylate free

Solvent free

Silicate free

EDTA free

PACKAGING:

Clothesline Fresh Iodine & Betadine Remover is packaged in RTU Handi Spray® quart bottles, 12 per case with 3 trigger sprayers. Label copy is available in English and Spanish. Secondary labels are also available.

Be sure to read all Directions, Precautionary and First Aid Statements on product labels before use of this or any Spartan product. If questions remain, consult your employer or a physician. Material Safety Data Sheets for all Spartan products are available from your authorized Spartan distributor or by visiting www.spartanchemical.com.

GUARANTEE:

Spartan's modern manufacturing and laboratory control insure uniform quality. If dissatisfied with performance of product, any unused portion may be returned for credit within one year of date of manufacture.